


Merlin HD

The Merlin HD is a very powerful digital magnifying glass. It useful for:

- **Low vision readers.**
- **Normal vision readers who need to read print that is small, faint, or otherwise hard to see.**
- **Anyone who wants to get a really close look at something.**

Instructions

- 1. Turn on the unit by pressing the Power switch .**
- 2. If the light does not turn on automatically, press the Light switch.**
- 3. Place what you want to view face-up on the platform.**
- 4. Turn the Size knob to the right to increase the level of magnification and to the left to decrease it.**
- 5. Slide the Brightness indicator to the right to increase the brightness and to the left to decrease it.**
- 6. Press the Mode button to change the appearance of the image. (If you change the Mode, it may also be necessary to adjust the Brightness.)**
- 7. Move the platform in any direction to scan the text. When the lever on the front of the platform is pushed to the right, the platform locks in place. To release the lock, push the lever to the right again.**

Note: Please ensure that the Merlin HD is turned off when it is not in use.